


PRODUCT CATALOGUE

Biomimetic **CORAL**




CORAL is the first implant system developed by Avinent and came about as a result of intensive research in the fields of biomedicine and materials engineering. CORAL implants have proven their outstanding predictability and supreme reliability during the years they have been present on the market. Their remarkable shape and innovative Biomimetic surface give excellent aesthetic results.

CONTENTS




Surgical phase	6
EC External connection	6
IC Internal connection	8
Accessories	10
Drill stops	12
Instruments and surgical kit	14
Drilling speed	15
Drill-bit length and marking	15
Avinent drilling protocol	16
Sterilization and packaging	18
Avinent Guide Surgery	20
Surgical motors	20
Radiographic template	20
Product index	21

Biomimetic CORAL



EXTERNAL CONNECTION

SURGICAL PHASE

	Platform Ø 3,5 mm					
	Ø 3,3 mm	Ref.	Ø 3,3 mm	Ref.	Ø 3,8 mm	Ref.
Implant 	-	-	-	-	3,8 x 7 mm	0836
	-	-	-	-	3,8 x 8,5 mm	0122
	3,3 x 10 mm	0123	3,3 x 10 mm	0184	3,8 x 10 mm	0121
	3,3 x 11,5 mm	0124	3,3 x 11,5 mm	0185	3,8 x 11,5 mm	0126
	3,3 x 13 mm	0125	3,3 x 13 mm	0186	3,8 x 13 mm	0127
	-	-	3,3 x 15 mm	0187	3,8 x 15 mm	0128
Healing abutment 	4 x 3 mm	0030	4 x 3 mm	3363	4 x 3 mm	3363
	4 x 4 mm	0031	4 x 5 mm	3364	4 x 5 mm	3364
	4 x 5 mm	0032	4 x 7 mm	3365	4 x 7 mm	3365
	5 x 3 mm	0033	5 x 3 mm	0036	5 x 3 mm	0036
	5 x 4 mm	0034	5 x 4 mm	0037	5 x 4 mm	0037
	5 x 5 mm	0035	5 x 5 mm	0003	5 x 5 mm	0003
	-	-	5 x 7 mm	0038	5 x 7 mm	0038
Anatomic healing abutment 	5 x 3 mm	0171	5 x 3 mm	0174	5 x 3 mm	0174
	5 x 4 mm	0172	5 x 4 mm	0175	5 x 4 mm	0175
	5 x 5 mm	0173	5 x 5 mm	0176	5 x 5 mm	0176
	-	-	5 x 7 mm	0177	5 x 7 mm	0177




Platform Ø 4,1 mm				Platform Ø 5,1 mm			
Ø 4,0 mm	Ref.	Ø 4,2 mm	Ref.	Ø 4,8 mm	Ref.	Ø 4,8 mm	Ref.
4,0 x 7 mm	0837	4,2 x 7 mm	0838	4,8 x 7 mm	0840	4,8 x 7 mm	0839
4,0 x 8,5 mm	0129	4,2 x 8,5 mm	0133	4,8 x 8,5 mm	0560	4,8 x 8,5 mm	0138
4,0 x 10 mm	0130	4,2 x 10 mm	0134	4,8 x 10 mm	0561	4,8 x 10 mm	0139
4,0 x 11,5 mm	0001	4,2 x 11,5 mm	0135	4,8 x 11,5 mm	0562	4,8 x 11,5 mm	0140
4,0 x 13 mm	0131	4,2 x 13 mm	0136	4,8 x 13 mm	0563	4,8 x 13 mm	0141
4,0 x 15 mm	0132	4,2 x 15 mm	0137	-	-	-	-
4 x 3 mm	3363	4 x 3 mm	3363	4 x 3 mm	3363	6 x 3 mm	0039
4 x 5 mm	3364	4 x 5 mm	3364	4 x 5 mm	3364	6 x 4 mm	0040
4 x 7 mm	3365	4 x 7 mm	3365	4 x 7 mm	3365	6 x 5 mm	0041
5 x 3 mm	0036	5 x 3 mm	0036	5 x 3 mm	0036	-	-
5 x 4 mm	0037	5 x 4 mm	0037	5 x 4 mm	0037	-	-
5 x 5 mm	0003	5 x 5 mm	0003	5 x 5 mm	0003	-	-
5 x 7 mm	0038	5 x 7 mm	0038	5 x 7 mm	0038	-	-
5 x 3 mm	0174	5 x 3 mm	0174	5 x 3 mm	0174	6 x 3 mm	0178
5 x 4 mm	0175	5 x 4 mm	0175	5 x 4 mm	0175	6 x 4 mm	0179
5 x 5 mm	0176	5 x 5 mm	0176	5 x 5 mm	0176	6 x 5 mm	0180
5 x 7 mm	0177	5 x 7 mm	0177	5 x 7 mm	0177	-	-

Biomimetic CORAL












INTERNAL CONNECTION












SURGICAL PHASE









		Platform Ø 3,5 mm			
		Ø 3,3 mm	Ref.	Ø 3,3 mm	Ref.
Implant 	-	-	-	-	
	-	-	-	-	
	3,3 x 10 mm	0385	3,3 x 10 mm	0408	
	3,3 x 11,5 mm	0386	3,3 x 11,5 mm	0409	
	3,3 x 13 mm	0387	3,3 x 13 mm	0410	
	3,3 x 15 mm	0388	3,3 x 15 mm	0411	
Healing abutment 	3,5 x 3 mm	3370	3,5 x 3 mm	3370	
	3,5 x 5 mm	3371	3,5 x 5 mm	3371	
	3,5 x 7 mm	3372	3,5 x 7 mm	3372	
	4,5 x 2 mm	0465	4,5 x 2 mm	0465	
	4,5 x 3 mm	0466	4,5 x 3 mm	0466	
	4,5 x 4 mm	0467	4,5 x 4 mm	0467	
	4,5 x 5 mm	0468	4,5 x 5 mm	0468	
	4,5 x 7 mm	0469	4,5 x 7 mm	0469	
	5,5 x 3 mm	0471	5,5 x 3 mm	0471	
	5,5 x 4 mm	0472	5,5 x 4 mm	0472	
	5,5 x 5 mm	0473	5,5 x 5 mm	0473	
	5,5 x 7 mm	0474	5,5 x 7 mm	0474	
	6,5 x 3 mm	2737	6,5 x 3 mm	2737	
	6,5 x 4 mm	2727	6,5 x 4 mm	2727	
6,5 x 5 mm	2726	6,5 x 5 mm	2726		
Anatomic healing abutment 	4,5 x 2 mm	0607	4,5 x 2 mm	0607	
	4,5 x 3 mm	0608	4,5 x 3 mm	0608	
	4,5 x 4 mm	0609	4,5 x 4 mm	0609	
	4,5 x 5 mm	0610	4,5 x 5 mm	0610	
	4,5 x 6 mm	0611	4,5 x 6 mm	0611	
	4,5 x 7 mm	0612	4,5 x 7 mm	0612	

Platform Ø 4,1 mm							
Ø 3,8 mm	Ref.	Ø 4,0 mm	Ref.	Ø 4,2 mm	Ref.	Ø 4,8 mm	Ref.
3,8 x 7 mm	0841	4,0 x 7 mm	0842	4,2 x 7 mm	0843	4,8 x 7 mm	0844
3,8 x 8,5 mm	0389	4,0 x 8,5 mm	0394	4,2 x 8,5 mm	0399	4,8 x 8,5 mm	0404
3,8 x 10 mm	0390	4,0 x 10 mm	0395	4,2 x 10 mm	0400	4,8 x 10 mm	0405
3,8 x 11,5 mm	0391	4,0 x 11,5 mm	0396	4,2 x 11,5 mm	0401	4,8 x 11,5 mm	0406
3,8 x 13 mm	0392	4,0 x 13 mm	0397	4,2 x 13 mm	0402	4,8 x 13 mm	0407
3,8 x 15 mm	0393	4,0 x 15 mm	0398	4,2 x 15 mm	0403	-	-
3,5 x 3 mm	3370	3,5 x 3 mm	3370	3,5 x 3 mm	3370	3,5 x 3 mm	3370
3,5 x 5 mm	3371	3,5 x 5 mm	3371	3,5 x 5 mm	3371	3,5 x 5 mm	3371
3,5 x 7 mm	3372	3,5 x 7 mm	3372	3,5 x 7 mm	3372	3,5 x 7 mm	3372
4,5 x 2 mm	0465	4,5 x 2 mm	0465	4,5 x 2 mm	0465	4,5 x 2 mm	0465
4,5 x 3 mm	0466	4,5 x 3 mm	0466	4,5 x 3 mm	0466	4,5 x 3 mm	0466
4,5 x 4 mm	0467	4,5 x 4 mm	0467	4,5 x 4 mm	0467	4,5 x 4 mm	0467
4,5 x 5 mm	0468	4,5 x 5 mm	0468	4,5 x 5 mm	0468	4,5 x 5 mm	0468
4,5 x 7 mm	0469	4,5 x 7 mm	0469	4,5 x 7 mm	0469	4,5 x 7 mm	0469
5,5 x 3 mm	0471	5,5 x 3 mm	0471	5,5 x 3 mm	0471	5,5 x 3 mm	0471
5,5 x 4 mm	0472	5,5 x 4 mm	0472	5,5 x 4 mm	0472	5,5 x 4 mm	0472
5,5 x 5 mm	0473	5,5 x 5 mm	0473	5,5 x 5 mm	0473	5,5 x 5 mm	0473
5,5 x 7 mm	0474	5,5 x 7 mm	0474	5,5 x 7 mm	0474	5,5 x 7 mm	0474
6,5 x 3 mm	2737	6,5 x 3 mm	2737	6,5 x 3 mm	2737	6,5 x 3 mm	2737
6,5 x 4 mm	2727	6,5 x 4 mm	2727	6,5 x 4 mm	2727	6,5 x 4 mm	2727
6,5 x 5 mm	2726	6,5 x 5 mm	2726	6,5 x 5 mm	2726	6,5 x 5 mm	2726
4,5 x 2 mm	0607	4,5 x 2 mm	0607	4,5 x 2 mm	0607	4,5 x 2 mm	0607
4,5 x 3 mm	0608	4,5 x 3 mm	0608	4,5 x 3 mm	0608	4,5 x 3 mm	0608
4,5 x 4 mm	0609	4,5 x 4 mm	0609	4,5 x 4 mm	0609	4,5 x 4 mm	0609
4,5 x 5 mm	0610	4,5 x 5 mm	0610	4,5 x 5 mm	0610	4,5 x 5 mm	0610
4,5 x 6 mm	0611	4,5 x 6 mm	0611	4,5 x 6 mm	0611	4,5 x 6 mm	0611
4,5 x 7 mm	0612	4,5 x 7 mm	0612	4,5 x 7 mm	0612	4,5 x 7 mm	0612

ACCESSORIES


		Ref.
Sterilization cassette OCEAN		1763
MiniBox		0526
Torque wrench		0295
Holding tool EC/IC		0651
Screwdrivers		Screwdriver ISO 1797 S (048) 0644
		Screwdriver ISO 1797 L (048) 0277
		Screwdriver ISO 1797 XL (048) 0645
		Implant driver ISO 1797 EC/IC S (2.5) 2693
		Implant driver ISO 1797 EC/IC L (2.5) 0300
		Implant driver ISO 1797 W&H EC/IC (2.5) 2692
		Screwdriver gold screw ISO 1797 0263
		Screwdriver Rhein 83 abutment
		Screwdriver grooved screw ISO 1797 0267
		Screwdriver impression coping closed tray ISO 1797 0723
		Screwdriver transepithelial abutment ISO 1797 S 0328
		Screwdriver transepithelial abutment ISO 1797 L 0726
		Screwdriver transepithelial angled abutment ISO 1797 S 0804
		Screwdriver transepithelial angled abutment ISO 1797 L 0648
		Implant driver ISO 1797 CC 3,5 S 2981
		Implant driver ISO 1797 CC 3,5 L 2984
		Implant driver ISO 1797 CC 4,1 S 2982
		Implant driver ISO 1797 CC 4,1 L 2985

Implant handle		Implant handle S Implant handle L		1878 2891	
Monoblock screwdrivers		Screwdriver S (048) Screwdriver L (048)		0274 0275	
		Implant driver S (2.5) Implant driver L (2.5)		0278 0299	
Drills		Pointed drill		0188	
			Pilot drill	2,0 - 2,8 mm	0094
				Twist drill	2,0 x 7 - 15 mm
		2,8 x 7 - 15 mm			0090
		3,2 x 7 - 15 mm			0091
		3,6 x 7 - 15 mm			0092
	Short twist drill	4,2 x 7 - 13 mm	0093		
		2,0 x 7 - 10 mm	0868		
		2,8 x 7 - 10 mm	0869		
		3,2 x 7 - 10 mm	0870		
	Hard bone twist drill	3,6 x 7 - 10 mm	0871		
		4,2 x 7 - 10 mm	0872		
		3,4 x 7 - 15 mm	0873		
		3,9 x 7 - 15 mm	0888		
Countersink			4,5 x 7 - 13 mm	0889	
			3,3 mm	0095	
			3,8 mm	0096	
			4,0 mm	0096	
			4,2 mm	0097	
Screw tap			4,8 mm	0098	
			3,3 mm	0099	
			3,8 mm	0100	
			4,0 mm	0100	
			4,2 mm	0101	
Direction indicator			4,8 mm	0102	
			1,9 - 2,7 mm	0021	
Drill extension ISO 1797				KI589B204	

Drill stops		Ø1,3-2,4 L7	6984
		Ø1,3-2,4 L8,5	6985
		Ø1,3-2,4 L10	6986
		Ø1,3-2,4 L11,5	6987
		Ø1,3-2,4 L13	6988
		Ø1,3-2,4 L15	6989
	Ø2,8-3,3 L7	6990	
	Ø2,8-3,3 L8,5	6991	
	Ø2,8-3,3 L10	6992	
	Ø2,8-3,3 L11,5	6993	
	Ø2,8-3,3 L13	6994	
	Ø2,8-3,3 L15	6995	
	Ø3,2-3,8 L7	6996	
	Ø3,2-3,8 L8,5	6997	
	Ø3,2-3,8 L10	6998	
	Ø3,2-3,8 L11,5	6999	
	Ø3,2-3,8 L13	7000	
	Ø3,2-3,8 L15	7001	
	Ø3,6-4,3 L7	7002	
	Ø3,6-4,3 L8,5	7003	
	Ø3,6-4,3 L10	7004	
	Ø3,6-4,3 L11,5	7005	
	Ø3,6-4,3 L13	7006	
	Ø3,6-4,3 L15	7007	
	Ø4,2-4,8 L7	7008	
	Ø4,2-4,8 L8,5	7009	
	Ø4,2-4,8 L10	7010	
	Ø4,2-4,8 L11,5	7011	
	Ø4,2-4,8 L13	7012	
	Ø5,7 L7	7013	
	Ø5,7 L8,5	7014	
	Ø5,7 L10	7015	
	Ø5,7 L11,5	7016	
Drill Stops Kit for implants up to Ø4,5		7959	
Drill Stops from Ø4,8 implants to Ø6,0		7960	


Instruments and surgical kit

Avinent supplies high-precision instruments designed specifically for inserting implants and suitable for use with internal and external connections. Avinent's two surgical boxes allow practitioners to choose the working system that best meets their needs.



The sterilization cassette CORAL stands out for its clear layout and attractive design, making it easy to use. The drilling sequence is clearly indicated by means of a simple colour code according to the diameter of the selected implant. All the pieces are placed in a sterilizable tray with a seethrough lid, giving a clear view of the interior.

The MiniBox is versatile, as it allows practitioners to select a specific sequence and take everything required for the surgical procedure with them in a small container. The box is sterilizable and can hold all the items needed for inserting implants.



Advisable torque

TYPE	VALUE	
Mechanical	35 Ncm	Screw for single/Multiple abutment*
	30 Ncm	Screw for angled titanium base**
	30 Ncm	Screw for angulation correction G2 (max. 30°)**
	20 Ncm	Screw for angulation correction G1 (max. 20°)**
	15 Ncm	Screw transepithelial angled abutment
	15 Ncm	Prosthetic screw for transepithelial
	35 Ncm	Transepithelial abutment (Uniblock / 2 parts)
	30 Ncm	LOCATOR® / LOCATOR R-Tx® Abutment
	25 Ncm	RHEIN83® Abutment (OT Equator)
Manual	≈8-15 Ncm	Scan Abutment
		Impression coping open/closed tray engaging
		Healing Abutment Ti
		Healing cap

* Regarding: Titanium base, Cemented abutment and Cemented angled abutment, Castable CoCr Base. Included: Gold screw.


** Except Transep. 4.8 of M1.40 which is 15 Ncm.

Drilling speed

	rpm
Guide drill	800 - 1.200
Drill ø 2 mm	800 - 1.200
Pilot drill	600 - 800
Drill ø 2,8 mm	400 - 600
Drill ø 3,2 mm	200 - 400
Drill ø 3,4 mm	200 - 400
Drill ø 3,6 mm	200 - 400
Drill ø 3,9 mm	200 - 400
Drill ø 4,2 mm	200 - 400
Drill ø 4,5 mm	200 - 400
Countersink	150
Screw Tap	20

Maximum recommended torque for implant insertion: 45-50 Ncm

Maximum recommended speed for implant insertion: 20 rpm



DRILL	A
ø 2,8 mm	0,70 mm
ø 3,2 mm	0,85 mm
ø 3,6 mm	0,85 mm
ø 4,2 mm	0,90 mm

Drill-bit length and marking

Avinent drills carry physical and laser markings to improve visibility during osteotomy, and follow a colour code according to the diameter of the implant. A two-part distinctive mark indicates those drill bits which are for hard bone. The marking corresponds to the length of the implant in crestal placement, but the distances are not absolute from the tip of the instrument to the mark. The length of the drill tip is not included in the depth mark, and so this distance must be taken into account when planning treatment and in carrying out the osteotomy.

* For any additional information and instrument maintenance instructions, please check www.avinent.com

AVINENT drilling protocol

The Avinent drilling protocol for inserting CORAL system implants is based on the anatomical and morphological characteristics of the bone to be worked on, the aim being to achieve satisfactory primary stability and a good bone-implant contact percentage, key factors for good osseointegration. The diameter of the helicoidal drills in the system is related to the core of the implant to be inserted.

Hence it is advisable to use the system's drilling sequence and instruments. The finish of the drills makes it easier to locate the marks that indicate length during surgery. Due to the fact that the surface treatment of the Avinent implant extends to the level of the platform, it is advisable to leave the implant at juxta-bone level.

* Hard bone drills are identified with two color indicators.

* The use of countersink is not indicated for 7 mm implants.

Implant \varnothing 3,3 mm

Ref.	0188	0089	0094	0090	0095	0099	
	BONE TYPE III-IV						
	BONE TYPE I-II						

Implant \varnothing 3,8/4,0 mm

Ref.	0188	0089	0094	0090	0091	0873	0096	0100	
	BONE TYPE III-IV						Option 1	Option 2	
	BONE TYPE I-II								

Implant ø 4,2 mm

Ref.	0188	0089	0094	0090	0091	0092	0888	0097	0101	
	BONE TYPE III-IV						Option 1	Option 2		
	BONE TYPE I-II									

Implant ø 4,8 mm

Ref.	0188	0089	0094	0090	0091	0092	0093	0889	0098	0102	
	BONE TYPE III-IV						Option 1	Option 2			
	BONE TYPE I-II										

Sterilization and packaging



EC

IC

Avinent implants are subjected to a sterilization process in accordance with the requirements of the CE mark for medical device. Avinent supplies its implants in packaging that is easily identified by a simple colour code.

The vial cap comes in different colours according to the diameter of the implant and has a sticker indicating the type of connection and the diameter and length of the implant.

The Avinent implant system is sold in a sterile blister pack, which ensures that the implant is fully protected until the given expiry date, so long as it is stored in suitable conditions.

Avinent supplies labeling with all its products that ensure every item can be properly traced by means of stickers that can be used in the patient's clinical history or in any other documentation required.



Blister pack (front and back)



Vial



Opening the vial



Opened vial

Opening the packaging

- Open the box and take out the blister pack.
- Open the blister pack and leave the vial containing the implant in the sterile working area, avoiding contact with any item that is not sterile. The sterile blister pack must not be opened prior to use.
- Open the vial by gripping the coloured cap with one hand while holding the transparent part in the other.
- Fit the implant driver (manual or hand piece) into the implant holder.
- Press to make a tight fit and raise the implant set slightly.
- Open the black cap on the bottom of the vial to remove the cover screw.

Information on the box

Rx only

CE 0197 **Biomimetic** CORAL

Manufacturer:
AVINENT Implant System S.L.U.
Pol. Ind. Santa Anna 1, 08251 Santpedor (Barcelona) Spain.

STERILE R



ifu.avinent.com



Implant box label

Implant EC
ø 3.3 x 10 (3.5)

REF 0123

ES Implante
PT Implante
FR Implant

(01) 08435452001046
(10) XXXXX
(17) AAMMDD

LOT XXXXX

use by AAAA-MM-DD



Labels that make
the product traceable

AVINENT
Implant EC ø 3.3 x 10 (3.5)

REF 0123 LOT XXXX

GTIN 08435452001046

STERILE R AAAA-MM-DD



AVINENT
Implant EC ø 3.3 x 10 (3.5)

REF 0123 LOT XXXX

GTIN 08435452001046

STERILE R AAAA-MM-DD



AVINENT
Implant EC ø 3.3 x 10 (3.5)

REF 0123 LOT XXXX

GTIN 08435452001046

STERILE R AAAA-MM-DD




AVINENT
Purchasing / Compra0123

REF 0123

LOT XXXXX

Implant EC/IC	External/internal connection indicator
ø 3.3 x 10 (3.5)	Diameter x length (implant holder)
REF 0123	Reference number
LOT XXXXX	Batch code
AAAA-MM-DD	Expiry date
STERILE R	Sterile. Sterilization method: radiation
ifu.avinent.com	Consult instructions for use
Manufacturer	Manufacturer
CE 0197	CE Mark NB num 0197. TÜV Rheinland
Rx only	(For the USA only) CAUTION: Federal law restricts this device to sale by or on the order of a physician
(01) 08435452001046 (10) XXXXX (17) AAMMDD	UDI (Unique Device Identifier Carrier)
Do not re-use	Do not re-use
Do not resterilize	Do not resterilize

Avinent Guided Surgery



Avinent Guided Surgery

- Surgical Guide Box 3DIEMME®
- Implant planning software with Real Guide®

CT or CBCT Scanners. Available libraries for:


- NewTom
- Carestream
- (More coming soon)



Surgical motors

Avinent offers a wide range of surgical motors designed for many clinical applications, produced by leading brands in the industry. Intended to simplify daily practice and for multiple users, they can be adapted to suit each practitioner's needs by allowing the entire drilling sequence to be customised.

Our surgical motors are ergonomically designed, easy to use and significantly help to reduce tiredness. They give implantology specialists a surgical unit for everyday use in safely carrying out oral procedures. All the surgical motors have been designed with the aim of minimising size and weight, to improve the balance of the instrument and motor in the practitioner's hand and thereby reduce tiredness during long treatments.



Radiographic template

The radiographic template is a guide to help practitioners precisely select the appropriate implant diameter and length. To adapt it to the proportions of the x-ray, the dimensions of the implant are shown at a scale 1:1 and 1.25:1.



Product index

REF.	DESCRIPTION	PAGE	REF.	DESCRIPTION	PAGE
0001	Implant EC 4,0 x 11,5 (4,1)	6	0139	Implant EC 4,8 x 10 (5,1)	6
0003	Healing abutment EC 4,1 x 5 x 5	6	0140	Implant EC 4,8 x 11,5 (5,1)	6
0021	Direction indicator ø 1,9 - ø 2,7	6	0141	Implant EC 4,8 x 13 (5,1)	6
0030	Healing abutment EC 3,5 x 4 x 3	6	0171	Anatomic healing abutment EC 3,5 x 5 x 3	6
0031	Healing abutment EC 3,5 x 4 x 4	6	0172	Anatomic healing abutment EC 3,5 x 5 x 4	6
0032	Healing abutment EC 3,5 x 4 x 5	6	0173	Anatomic healing abutment EC 3,5 x 5 x 5	6
0033	Healing abutment EC 3,5 x 5 x 3	6	0174	Anatomic healing abutment EC 4,1 x 5 x 3	6
0034	Healing abutment EC 3,5 x 5 x 4	6	0175	Anatomic healing abutment EC 4,1 x 5 x 4	6
0035	Healing abutment EC 3,5 x 5 x 5	6	0176	Anatomic healing abutment EC 4,1 x 5 x 5	6
0036	Healing abutment EC 4,1 x 5 x 3	6	0177	Anatomic healing abutment EC 4,1 x 5 x 7	6
0037	Healing abutment EC 4,1 x 5 x 4	6	0178	Anatomic healing abutment EC 5,1 x 6 x 3	6
0038	Healing abutment EC 4,1 x 5 x 7	6	0179	Anatomic healing abutment EC 5,1 x 6 x 4	6
0039	Healing abutment EC 5,1 x 6 x 3	6	0180	Anatomic healing abutment EC 5,1 x 6 x 5	6
0040	Healing abutment EC 5,1 x 6 x 4	6	0184	Implant EC 3,3 x 10 (4,1)	6
0041	Healing abutment EC 5,1 x 6 x 5	6	0185	Implant EC 3,3 x 11,5 (4,1)	6
0082	Screwdriver ISO 1797 (035)	10	0186	Implant EC 3,3 x 13 (4,1)	6
0089	Twist drill ø 2,0 x 7,0 - 15	11	0187	Implant EC 3,3 x 15 (4,1)	6
0090	Twist drill ø 2,8 x 7,0 - 15	11	0188	Pointed drill	11
0091	Twist drill ø 3,2 x 7,0 - 15	11	0263	Screwdriver gold screw ISO 1797 Screwdriver Rhein 83 abutment	10
0092	Twist drill ø 3,6 x 7,0 - 15	11	0267	Screwdriver grooved screw ISO 1797	10
0093	Twist drill ø 4,2 x 7,0 - 13	11	0274	Screwdriver S (048)	10
0094	Pilot drill ø 2,0 - ø 2,8 x 7	11	0275	Screwdriver L (048)	10
0095	Countersink 3,3	11	0277	Screwdriver ISO 1797 L (048)	10
0096	Countersink 3,8 - 4,0	11	0278	Implant driver S (2,5)	10
0097	Countersink 4,2	11	0295	Torque wrench	10
0098	Countersink 4,8	11	0296	Sterilization cassette CORAL	10
0099	Screw tap 3,3	11	0299	Implant driver L (2,5)	10
0100	Screw tap 3,8 - 4,0	11	0300	Implant driver ISO 1797 EC/IC L (2,5)	10
0101	Screw tap 4,2	11	0328	Screwdriver transepithelial abutment ISO 1797 S	10
0102	Screw tap 4,8	11	0385	Implant IC 3,3 x 10 (3,5)	8
0121	Implant EC 3,8 x 10 (4,1)	6	0386	Implant IC 3,3 x 11,5 (3,5)	8
0122	Implant EC 3,8 x 8,5 (4,1)	6	0387	Implant IC 3,3 x 13 (3,5)	8
0123	Implant EC 3,3 x 10 (3,5)	6	0388	Implant IC 3,3 x 15 (3,5)	8
0124	Implant EC 3,3 x 11,5 (3,5)	6	0389	Implant IC 3,8 x 8,5 (4,1)	8
0125	Implant EC 3,3 x 13 (3,5)	6	0390	Implant IC 3,8 x 10 (4,1)	8
0126	Implant EC 3,8 x 11,5 (4,1)	6	0391	Implant IC 3,8 x 11,5 (4,1)	8
0127	Implant EC 3,8 x 13 (4,1)	6	0392	Implant IC 3,8 x 13 (4,1)	8
0128	Implant EC 3,8 x 15 (4,1)	6	0393	Implant IC 3,8 x 15 (4,1)	8
0129	Implant EC 4,0 x 8,5 (4,1)	6	0394	Implant IC 4,0 x 8,5 (4,1)	8
0130	Implant EC 4,0 x 10 (4,1)	6	0395	Implant IC 4,0 x 10 (4,1)	8
0131	Implant EC 4,0 x 13 (4,1)	6	0396	Implant IC 4,0 x 11,5 (4,1)	8
0132	Implant EC 4,0 x 15 (4,1)	6	0397	Implant IC 4,0 x 13 (4,1)	8
0133	Implant EC 4,2 x 8,5 (4,1)	6	0398	Implant IC 4,0 x 15 (4,1)	8
0134	Implant EC 4,2 x 10 (4,1)	6	0399	Implant IC 4,2 x 8,5 (4,1)	8
0135	Implant EC 4,2 x 11,5 (4,1)	6	0400	Implant IC 4,2 x 10 (4,1)	8
0136	Implant EC 4,2 x 13 (4,1)	6	0401	Implant IC 4,2 x 11,5 (4,1)	8
0137	Implant EC 4,2 x 15 (4,1)	6			
0138	Implant EC 4,8 x 8,5 (5,1)	6			

REF.	DESCRIPTION	PAGE	REF.	DESCRIPTION	PAGE
0402	Implant IC 4,2 x 13 (4,1)	8	0839	Implant EC 4,8 x 7 (5,1)	6
0403	Implant IC 4,2 x 15 (4,1)	8	0840	Implant EC 4,8 x 7 (4,1)	6
0404	Implant IC 4,8 x 8,5 (4,1)	8	0841	Implant IC 3,8 x 7 (4,1)	8
0405	Implant IC 4,8 x 10 (4,1)	8	0842	Implant IC 4,0 x 7 (4,1)	8
0406	Implant IC 4,8 x 11,5 (4,1)	8	0843	Implant IC 4,2 x 7 (4,1)	8
0407	Implant IC 4,8 x 13 (4,1)	8	0844	Implant IC 4,8 x 7 (4,1)	8
0408	Implant IC 3,3 x 10 (4,1)	8	0868	Twist drill ø 2,0 x 7,0 - 10	11
0409	Implant IC 3,3 x 11,5 (4,1)	8	0869	Twist drill ø 2,8 x 7,0 - 10	11
0410	Implant IC 3,3 x 13 (4,1)	8	0870	Twist drill ø 3,2 x 7,0 - 10	11
0411	Implant IC 3,3 x 15 (4,1)	8	0871	Twist drill ø 3,6 x 7,0 - 10	11
0465	Anatomic healing abutment IC 3,5/4,1 x 4,5 x 2	8	0872	Twist drill ø 4,2 x 7,0 - 10	11
0466	Anatomic healing abutment IC 3,5/4,1 x 4,5 x 3	8	0873	Twist drill ø 3,4 x 7,0 - 15	11
0467	Anatomic healing abutment IC 3,5/4,1 x 4,5 x 4	8	0888	Twist drill ø 3,9 x 7,0 - 15	11
0468	Anatomic healing abutment IC 3,5/4,1 x 4,5 x 5	8	0889	Twist drill ø 4,5 x 7,0 - 15	11
0469	Anatomic healing abutment IC 3,5/4,1 x 4,5 x 7	8	1878	Implant handle S	10
0471	Anatomic healing abutment IC 3,5/4,1 x 5,5 x 3	8	2692	Implant driver ISO 1797 W&H EC/IC (2,5)	10
0472	Anatomic healing abutment IC 3,5/4,1 x 5,5 x 4	8	2693	Implant driver ISO 1797 EC/IC S (2,5)	10
0473	Anatomic healing abutment IC 3,5/4,1 x 5,5 x 5	8	2726	Anatomic healing abutment IC 3,5/4,1 x 6,5 x 5	8
0474	Anatomic healing abutment IC 3,5/4,1 x 5,5 x 7	8	2727	Anatomic healing abutment IC 3,5/4,1 x 6,5 x 4	8
0526	MiniBox	10	2737	Anatomic healing abutment IC 3,5/4,1 x 6,5 x 3	8
0527	Radel plate 3x6"	10	2891	Implant handle L	10
0560	Implant EC 4,8 x 8,5 (4,1)	8	3363	Healing abutment EC 4,1 x 4 x 3	6
0561	Implant EC 4,8 x 10 (4,1)	8	3364	Healing abutment EC 4,1 x 4 x 5	6
0562	Implant EC 4,8 x 11,5 (4,1)	8	3365	Healing abutment EC 4,1 x 4 x 7	6
0563	Implant EC 4,8 x 13 (4,1)	8	3366	Anatomic healing abutment PEEK EC 4,1 x 6 x 5	6
0607	Aesthetic healing abutment IC 4,5 x 2	8	3367	Anatomic healing abutment PEEK EC 4,1 x 8 x 5	6
0608	Aesthetic healing abutment IC 4,5 x 3	8	3368	Anatomic healing abutment PEEK EC 5,1 x 6 x 5	6
0609	Aesthetic healing abutment IC 4,5 x 4	8	3369	Anatomic healing abutment PEEK EC 5,1 x 8 x 5	6
0610	Aesthetic healing abutment IC 4,5 x 5	8	3370	Anatomic healing abutment IC 3,5/4,1 x 3,5 x 3	8
0611	Aesthetic healing abutment IC 4,5 x 6	8	3371	Anatomic healing abutment IC 3,5/4,1 x 3,5 x 5	8
0612	Aesthetic healing abutment IC 4,5 x 7	8	3372	Anatomic healing abutment IC 3,5/4,1 x 3,5 x 7	8
0644	Screwdriver S ISO 1797 (048)	10	6984	Drill Stop Ø1.3-2.4 L7	12
0645	Screwdriver XL ISO 1797 (048)	10	6985	Drill Stop Ø1.3-2.4 L8.5	12
0648	Screwdriver transepithelial angled abutment ISO 1797 L	10	6986	Drill Stop Ø1.3-2.4 L10	12
0651	Holding tool EC/IC	10	6987	Drill Stop Ø1.3-2.4 L11.5	12
KI589B204	Drill extension ISO 1797	10	6988	Drill Stop Ø1.3-2.4 L13	12
0723	Screwdriver coping closed tray ISO 1797	10	6989	Drill Stop Ø1.3-2.4 L15	12
0726	Screwdriver transepithelial abutment ISO 1797 L	10	6990	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L7	12
0790	Handle L ISO 1797	10	6991	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L8.5	12
0791	Handle S ISO 1797	10	6992	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L10	12
0804	Screwdriver transepithelial angled abutment ISO 1797 S	10	6993	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L11.5	12
0836	Implant EC 3,8 x 7 (4,1)	6	6994	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L13	12
0837	Implant EC 4,0 x 7 (4,1)	6	6995	CORAL/OCEAN/ICEBERG Drill Stop Ø2.8-3.3 L15	12
0838	Implant EC 4,2 x 7 (4,1)	6	6996	CORAL/OCEAN/ICEBERG Drill Stop Ø3.2-3.8 L7	12

*The availability of the products in this catalog may vary depending on the country where you are based. Please contact your Avinent distributor for further information.

REF.	DESCRIPTION	PAGE
7001	CORAL/OCEAN/ICEBERG Drill Stop Ø3.2-3.8 L15	12
7002	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L7	12
7003	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L8.5	12
7004	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L10	12
7005	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L11.5	12
7006	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L13	12
7007	CORAL/OCEAN/ICEBERG Drill Stop Ø3.6-4.3 L15	12
7008	CORAL/OCEAN/ICEBERG Drill Stop Ø4.2-4.8 L7	12
7009	CORAL/OCEAN/ICEBERG Drill Stop Ø4.2-4.8 L8.5	12
7010	CORAL/OCEAN/ICEBERG Drill Stop Ø4.2-4.8 L10	12
7011	CORAL/OCEAN/ICEBERG Drill Stop Ø4.2-4.8 L11.5	12
7012	CORAL/OCEAN/ICEBERG Drill Stop Ø4.2-4.8 L13	12
7013	OCEAN Drill Stop Ø5.7 L7	12
7014	OCEAN Drill Stop Ø5.7 L8.5	12
7015	OCEAN Drill Stop Ø5.7 L10	12
7016	OCEAN Drill Stop Ø5.7 L11.5	12
7959	Drill Stop Kit for implants up to Ø4.5	12
7960	Drill Stop Set for implants from Ø4.8 to Ø6.0	12





Carretera de Navarcles, 107
Pol. Industrial Santa Anna I
08251 Santpedor (Barcelona) - España
T. (+34) 93 827 34 65
www.avinent.com - avinent@avinent.com

AVINENT AUSTRALIA
australia@avinent.com

AVINENT BENELUX
benelux@avinent.com

AVINENT CANADA
canada@avinent.com

AVINENT COLOMBIA
colombia@avinent.com

AVINENT CROATIA
croatia@avinent.com

AVINENT FRANCE
france@avinent.com

AVINENT KUWAIT
kuwait@avinent.com

AVINENT LITHUANIA
lithuania@avinent.com

AVINENT POLAND
poland@avinent.com

AVINENT PORTUGAL
portugal@avinent.com

AVINENT ROMANIA
romania@avinent.com

AVINENT TAIWAN
taiwan@avinent.com

AVINENT UKRAINE
ukraine@avinent.com

AVINENT UNITED ARAB EMIRATES
uae@avinent.com

AVINENT UNITED KINGDOM
uk@avinent.com

AVINENT USA
usa@avinent.com

AVINENT VIETNAM
vietnam@avinent.com

AVINENT SAUDI ARABIA
ksa@avinent.com



Management
System
ISO 13485:2016
www.tuv.com

