

PRODUCT CATALOGUE

Biomimetic **PEARL**






PEARL is the Avinent system of mini implants, an innovative product with unique features that aims to provide a quality alternative to what is available today as regards the range of solutions for removable, minimally invasive prostheses.

CONTENTS

Biomimetic Pearl	4
Characteristics and specifications	6
Surgical phase	8
Prosthetic phase	8
Accessories	9
Drills stops	10
Drill-bit length and marking	12
Instruments and surgical kit	12
Drilling speed	12
AVINENT drilling protocol	13
Sterilization and packaging	14
Product index	16



Biomimetic PEARL


A unique and resistant mini implant

The Biomimetic PEARL system of mini implants is Avinent's innovative solution to provide a quality product as regards the range of systems for removable, minimally invasive prostheses.

It is a unique mini implant whose exclusive two-piece design makes it very easy to position the attachment and thereby control soft tissue. In addition, its polished collar is an ideal feature for unbeatable aesthetic results.

PEARL has been designed to help surgeons cover various gingival heights for a much more complete solution. It also benefits from a single protocol and unprecedented versatility, in addition to a very intuitive surgical procedure.

The PEARL mini implant also preserves the unique characteristics of CORAL and OCEAN implants as regards the BIOMIMETIC ADVANCED SURFACE (BAS), which is based on the biochemical processes taking place in nature that encourage osseointegration.



- ▶ Innovative two-piece design to help insert the attachment.
- ▶ Three different diameters, 2,0 / 2,4 / 2,8 mm, for a wide range of possibilities in all cases and maintaining the same platform in all diameters (3,0 mm).
- ▶ Polished collar for improved aesthetic results.
- ▶ Continuous progressive double-thread that facilitates implant insertion.
- ▶ BIOMIMETIC ADVANCED SURFACE, which increases bone-implant contact (BIC) and accelerates osseointegration.
- ▶ Easy and intuitive drilling protocol.



Unique two-piece design

PEARL is characterized by an innovative two-piece design (implant and abutment) that provides great versatility to practitioners and represents an ideal solution for removable, minimally invasive prostheses, as well as allowing immediate loading in cases of complete restoration and osseointegration over time.

The mini implant also features a polished collar for improved aesthetic results.

ø 2,0 mm



ø 2,4 mm



ø 2,8 mm




Three diameters for all solutions

The PEARL system has three different diameters providing a wide range of possibilities in all cases and maintaining the 3,0 mm platform in all diameters.


PEARL's 2,0 mm diameter is suitable for mandibles, its 2,4 mm for mandibles and maxillae and its 2,8 mm is recommended for maxillae, a perfect choice for all types of cases.

Intuitive, effective placement


The system incorporates a transporter whose function is to prevent soft tissue collapse from insertion of the implant to the placement of the abutment.




1
Placement of the implant alongside the transporter




2
Removal of the screw once the transporter is positioned




3
Placement of guided Rhein abutment thanks to the transporter, avoiding soft tissue collapse




4
Removal of transporter



5
Tissue fit around the abutment



6
Placement of final restoration



Attachment to control soft tissues

The PEARL system has three different attachments to resolve height differences in gums of varying thickness. This fact implies a safety margin for practitioners when it comes to positioning the mini implant because it provides them with different options according to the final placement result.

Biomimetic PEARL

SURGICAL PHASE

		Platform Ø 3 mm					
		Ø 2,0 mm	Ref.	Ø 2,4 mm	Ref.	Ø 2,8 mm	Ref.
Implant		2,0 x 10 mm	3819	2,4 x 10 mm	3917	2,8 x 10 mm	3920
		2,0 x 11,5 mm	3916	2,4 x 11,5 mm	3918	2,8 x 11,5 mm	3921
		2,0 x 13 mm	3826	2,4 x 13 mm	3919	2,8 x 13 mm	3922
		2,0 x 15 mm	3930	2,4 x 15 mm	3931	2,8 x 15 mm	3932
Cover screw		4199					
		Platform Ø 3 mm					
		Ø 3,5 mm	Ref.	Ø 4,0 mm	Ref.		
Healing abutment		H 2,5	6016	H 2,5	6019		
		H 3,5	6017	H 3,5	6020		
		H 5,5	6018	H 5,5	6021		

PROSTHETIC PHASE

		Platform Ø 3 mm	
		Ø 2,5 mm	Ref.
Overdenture abutment		2,5 x 0,8 mm 2,5 x 1,8 mm 2,5 x 2,8 mm	4166 4167 4168
Processing cap (2 u.) ¹		3317	
Processing cap Smart Box ²		4210	
Smart Box Housing ²		4211	
Tool for insertion and extraction caps		4209	
Impression coping		3318	
Analog		3319	
Retention replacement male purple** (2,7 kgf)		3322	
Retention replacement male white** (1,8 kgf)		3323	
Retention replacement male pink** (1,2 kgf)		3324	
Retention replacement male yellow** (0,6 kgf)		3325	
Retention replacement male black** (Only for laboratories)		4695	
Retention replacement Smart Box** (Only for laboratories)		4212	

ACCESSORIES


		Ref.
MiniBox		0526
Torque wrench		0295
Screwdrivers		Screwdriver ISO 1797 L (048) 0277
		Implant driver ISO 1797 EC/IC L (2.5) 2693
		Screwdriver gold screw ISO 1797 Screwdriver Rhein 83 abutment 0263
Screwdriver handle		Handle ISO 1797 S 0791
Implant handle		Handle ISO 1797 S Handle ISO 1797 L 1878 2891
Monoblock screwdrivers		Screwdriver S (048) 0274
		Implant driver S (2.5) 0278
Drills		Tissue punch 4165
		Drill 1,3 x 10 - 15 mm 3956 1,6 x 7 - 15 mm 2046 2,0 x 7 - 15 mm 0089 2,4 x 10 - 15 mm 3957
Drill extension ISO 1797		KI589B204

** In blisters of four units

¹ Maximum angulation 25°

² Maximum angulation 50°

Drill stops		Ø1,3-2,4 L7	6984
		Ø1,3-2,4 L8,5	6985
		Ø1,3-2,4 L10	6986
		Ø1,3-2,4 L11,5	6987
		Ø1,3-2,4 L13	6988
		Ø1,3-2,4 L15	6989
		Ø2,8-3,3 L7	6990
		Ø2,8-3,3 L8,5	6991
		Ø2,8-3,3 L10	6992
		Ø2,8-3,3 L11,5	6993
		Ø2,8-3,3 L13	6994
		Ø2,8-3,3 L15	6995
		Ø3,2-3,8 L7	6996
		Ø3,2-3,8 L8,5	6997
		Ø3,2-3,8 L10	6998
		Ø3,2-3,8 L11,5	6999
		Ø3,2-3,8 L13	7000
		Ø3,2-3,8 L15	7001
		Ø3,6-4,3 L7	7002
		Ø3,6-4,3 L8,5	7003
		Ø3,6-4,3 L10	7004
		Ø3,6-4,3 L11,5	7005
		Ø3,6-4,3 L13	7006
		Ø3,6-4,3 L15	7007
Drill Stops Kit for implants up to Ø4,5			7959



Drill-bit length and marking

Avinent drills carry laser markings to improve visibility during osteotomy and follow a point code according to the drilling sequence. The marking corresponds to the length of the implant in crestal placement, but the distances are not absolute from the tip of the instrument to the mark. The length of the drill tip is not included in the depth mark, so this distance must be taken into account when planning treatment and in carrying out the osteotomy.

DRILL	A
ø 1,3 mm	0,35 mm
ø 1,6 mm	0,35 mm
ø 2,0 mm	0,50 mm
ø 2,4 mm	0,50 mm


Advisable torque

TYPE	VALUE	
Mechanical	35 Ncm	Screw for single abutment*
	15 Ncm	Temporary Abutment Ti
	25 Ncm	RHEIN83° Abutment (OT Equator)
Manual	≈8-15 Ncm	Scan Abutment
		Impression coping open/closed tray engaging
		Healing Abutment Ti
		Healing cap

* Regarding: Titanium base, Cemented abutment and Cemented angled abutment, Castable CoCr Base.

Instruments and surgical kit

Avinent provides specialists with high precision instruments whose design is adapted to the implant. Avinent's smaller and more versatile compact kit allows you to select a specific sequence and take everything required for the surgical procedure with you in a small container. The kit is sterilizable and can hold all the items needed for inserting prostheses.



Drilling speed

	rpm
Drill ø 1,3 mm	800 - 1.200
Drill ø 1,6 mm	200 - 400 *
Drill ø 2,0 mm	200 - 400 *
Drill ø 2,4 mm	200 - 400

* Important note: A drilling speed of 800 - 1.200 rpm is recommended when using 1,6 mm and 2,0 mm drill bits for initial drilling.


Maximum recommended torque for implant insertion: 45 Ncm

Recommended torque to achieve a good primary stability: 30-35 Ncm

Maximum recommended speed for implant insertion: 20 rpm

Tissue punch

The Avinent tool that can make a circular cut of soft tissue. The markings it carries help to identify which Rhein abutment needs to be placed according to gum thickness.



Mark 1*, 2* and 3 mm:
suitable for 0,8 mm attachment

Mark 4 mm:
suitable for 1,8 mm attachment


Mark 5 mm:
suitable for 2,8 mm attachment


* Marks 1 and 2 mm will leave the implant exposed.


Avinent drilling protocol

The PEARL system surgical drilling protocol makes for easy, intuitive and streamlined placement in order to prioritize soft tissue and emergency profile preservation for excellent, practical results.

Avinent PEARL is suitable for all bone types. The diameter of the system's helical drill bits is related to the implant core to be inserted, so it is advisable to follow the system's own drilling and instrument sequence. The finish of the drills makes it easier to locate the length indicator marks during surgery.

				IMPLANT ø 2,0 mm	
					
Tissue punch ø 2,4 mm		Drill ø 1,3 mm	Drill ø 1,6 mm	Implant	
Ref. 4165		3956	2046		
BONE TYPE III-IV					
BONE TYPE I-II					

				IMPLANT ø 2,4 mm	
					
Tissue punch ø 2,4 mm		Drill ø 1,6 mm	Drill ø 2 mm	Implant	
Ref. 4165		2046	0089		
BONE TYPE III-IV					
BONE TYPE I-II					

				IMPLANT ø 2,8 mm	
					
Tissue punch ø 2,4 mm		Drill ø 2 mm	Drill ø 2,4 mm	Implant	
Ref. 4165		0089	3957		
BONE TYPE III-IV					
BONE TYPE I-II					

Sterilization and packaging



3.0



Avinent implants are subjected to a sterilization process in accordance with the requirements of the CE mark for medical products. Avinent supplies its implants in packaging that is easily identified by a simple color code. The vial stopper is white and has a sticker indicating the diameter and length of the implant.

The Avinent implant system is sold in a sterile blister pack, which ensures that the implant is fully protected until the given expiry date, provided that it is stored in suitable conditions.

Avinent supplies labeling with all its products that ensures every item can be properly traced by means of stickers that can be used in the patient's case history or in any other documentation require.




Blister pack (front and back)



Vial



Opening the vial



Implant

Opening the packaging

- Open the box and take out the blister pack.
- Open the blister pack and leave the vial containing the implant in the sterile working area, avoiding contact with any item that is not sterile. The sterile blister pack must not be opened prior to use.
- Open the vial by gripping the white stopper with one hand while holding the transparent part in the other.
- Fit the transporting tool (manual or mechanical for contra-angle) to the implant holder.
- Press to make a tight fit and raise the implant set slightly.

Information on implant box

Rx only

CE **Biomimetic** PEARL

0197
Manufacturer:
AVINENT Implant System S.L.U.
Pol. Ind. Santa Anna I, 08251 Santpedor (Barcelona) Spain.

STERILE R



ifu.avinent.



Implant
box label

Implant
ø 2.8 x 10 (3.0) Ti
REF 1558

ES Implant
PT Implant
FR Implant

(01) 08435452001046
(10) XXXXX
(17) AAMMDD

LOT XXXX

USE BY AAAA-MM-DD

Labels for product
traceability

AVINENT
Implant ø 2.8 x 10 (3.0)
REF 1558 LOT XXXX
GTIN 08435452001046
STERILE R AAAA-MM-DD

AVINENT
Implant ø 2.8 x 10 (3.0)
REF 1558 LOT XXXX
GTIN 08435452001046
STERILE R AAAA-MM-DD

AVINENT
Implant ø 2.8 x 10 (3.0)
REF 1558 LOT XXXX
GTIN 08435452001046
STERILE R AAAA-MM-DD

AVINENT
Purchasing / Compra
REF 1558
LOT XXXX

ø 2.8 x 10 (3.0)	Diameter x length (implant holder)
Ti	Titanium
REF 1558	Reference number
LOT xxxxx	Batch code
AAAA-MM-DD	Expiry date
STERILE R	Sterile. Sterilization method: radiation
ifu.avinent.com	Consult instructions for use
CE	Manufacturer
Rx only	CE Mark NB num 0197. TÜV Rheinland
(01) 08435452001046 (10) XXXXX (17) AAMMDD	(For the USA only) CAUTION: Federal law restricts this device to sale by or on the order of a physician
UDI	UDI (Unique Device Identifier Carrier)
Do not re-use	Do not re-use
Do not re-sterilize	Do not re-sterilize

Product index

REF.	DESCRIPTION	PAGE	REF.	DESCRIPTION	PAGE
0089	Drill ø 2,0 x 7 - 15 mm	W9	4695	Retention replacement male black (4 u.) (Only for laboratories)	8
0263	Screwdriver gold screw ISO 1797 Screwdriver Rhein 83 abutment	9	6016	Healing abutment H 2,5 ø 3,5 (3,0)	9
0274	Screwdriver S (048)	9	6017	Healing abutment H 3,5 ø 3,5 (3,0)	9
0277	Screwdriver ISO 1797 L (048)	9	6018	Healing abutment H 5,5 ø 3,5 (3,0)	9
0278	Implant driver S (2.5)	9	6019	Healing abutment H 2,5 ø 4,0 (3,0)	9
0295	Torque wrench	9	6020	Healing abutment H 3,5 ø 4,0 (3,0)	9
0526	MiniBox	9	6021	Healing abutment H 5,5 ø 4,0 (3,0)	9
K1589B204	Drill extension	9	6984	Drill Stop Ø1.3-2.4 L7	12
0791	Handle ISO 1797 S	9	6985	Drill Stop Ø1.3-2.4 L8.5	12
1878	Implant handle S	9	6986	Drill Stop Ø1.3-2.4 L10	12
2046	Drill ø 1,6 x 7 - 15 mm	9	6987	Drill Stop Ø1.3-2.4 L11.5	12
2693	Implant driver ISO 1797 EC/IC L (2.5)	9	6988	Drill Stop Ø1.3-2.4 L13	12
2891	Implant handle L	9	6989	Drill Stop Ø1.3-2.4 L15	12
3317	Rhein 83 Processing cap (2 u.)1	8	6990	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L7	12
3318	Rhein 83 impression coping	8	6991	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L8,5	12
3319	Rhein 83 analog (2 units)	8	6992	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L10	12
3322	Retention replacement male purple (2,7 kgf) (4 u.)	9	6993	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L11,5	12
3323	Retention replacement male transparent (1,8 kgf) (4 u.)	9	6994	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L13	12
3324	Retention replacement male pink (1,2 kgf) (4 u.)	9	6995	CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L15	12
3325	Retention replacement male yellow (0,6 kgf) (4 u.)	9	6996	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L7	12
3819	Biomimetic Pearl implant 2,0 x 10 (3,0)	8	6997	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L8,5	12
3826	Biomimetic Pearl implant 2,0 x 13 (3,0)	8	6998	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L10	12
3916	Biomimetic Pearl implant 2,0 x 11,5 (3,0)	8	6999	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L11,5	12
3917	Biomimetic Pearl implant 2,4 x 10 (3,0)	8	7000	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L13	12
3918	Biomimetic Pearl implant 2,4 x 11,5 (3,0)	8	7001	CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L15	12
3919	Biomimetic Pearl implant 2,4 x 13 (3,0)	8	7002	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L7	12
3920	Biomimetic Pearl implant 2,8 x 10 (3,0)	8	7003	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L8,5	12
3921	Biomimetic Pearl implant 2,8 x 11,5 (3,0)	8	7004	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L10	12
3922	Biomimetic Pearl implant 2,8 x 13 (3,0)	8	7005	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L11,5	12
3930	Biomimetic Pearl implant 2,0 x 15 (3,0)	8	7006	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L13	12
3931	Biomimetic Pearl implant 2,4 x 15 (3,0)	8	7007	CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L15	12
3932	Biomimetic Pearl implant 2,8 x 15 (3,0)	8	7959	Drill Stop Kit for implants up to Ø4.5	12
3956	Drill ø 1,3 x 10 - 15 mm	9			
3957	Drill ø 2,4 x 10 - 15 mm	9			
4165	Tissue punch	9			
4166	Rhein 83 abutment 2,5 x 0,8	8			
4167	Rhein 83 abutment 2,5 x 1,8	8			
4168	Rhein 83 abutment 2,5 x 2,8	8			
4199	Cover screw	8			
4209	Rhein 83 tool for insertion and extraction caps	8			
4210	Processing cap Smart Box	8			
4211	Smart Box Housing	8			
4212	Retention replacement Smart Box (4 u.) (Only for laboratories)	8			



Carretera de Navarcles, 107
Pol. Industrial Santa Anna I
08251 Santpedor (Barcelona) - España
T. (+34) 93 827 34 65
www.avinent.com - avinent@avinent.com

AVINENT AUSTRALIA
australia@avinent.com

AVINENT BENELUX
benelux@avinent.com

AVINENT CANADA
canada@avinent.com

AVINENT COLOMBIA
colombia@avinent.com

AVINENT CROATIA
croatia@avinent.com

AVINENT FRANCE
france@avinent.com

AVINENT KUWAIT
kuwait@avinent.com

AVINENT LITHUANIA
lithuania@avinent.com

AVINENT POLAND
poland@avinent.com

AVINENT PORTUGAL
portugal@avinent.com

AVINENT ROMANIA
romania@avinent.com

AVINENT TAIWAN
taiwan@avinent.com

AVINENT UKRAINE
ukraine@avinent.com

AVINENT UNITED ARAB EMIRATES
uae@avinent.com

AVINENT UNITED KINGDOM
uk@avinent.com

AVINENT USA
usa@avinent.com

AVINENT VIETNAM
vietnam@avinent.com

AVINENT SAUDI ARABIA
ksa@avinent.com



Management
System
ISO 13485:2016
www.tuv.com

