


PRODUCT CATALOGUE

Biomimetic **ICEBERG**

The **ICEBERG** system is a new line that adds to the range of Avinent implants and that is especially aimed at simplifying post-surgery treatments and increasing the creation of soft tissue around the implant thanks to its polished neck at a tissue level. The ICEBERG system offers a very simple and logical surgical procedure, with its own surgical kit and a simple and practical drilling protocol thanks to its color code.



CONTENTS




| | |
|--------------------------------------|----|
| Surgical phase | 4 |
| CCI Iceberg Connection | 4 |
| Accessories | 6 |
| Drill Stops | 8 |
| Instruments and surgical kit | 10 |
| Drilling speed | 11 |
| Drill-bit length and marking | 11 |
| Accessories CC.I | 11 |
| Platform indicators | 11 |
| Avinent drilling biological protocol | 12 |
| Avinent drilling standard protocol | 14 |
| Sterilization and packaging | 16 |
| Avinent Guided Surgery | 18 |
| Micromotors | 18 |
| Radiographic template | 18 |
| Product index | 19 |

Biomimetic ICEBERG



























ICEBERG CONNECTION



SURGICAL PHASE









| | Platform Ø 3,5 mm | | | | | | Platform Ø 4,1 mm | | | | | |
|---|-------------------|------|---------------|------|---------------|------|-------------------|------|---------------|------|---------------|------|
| | Ø 3,5 mm | Ref. | Ø 4,0 mm | Ref. | Ø 4,5 mm | Ref. | Ø 5,0 mm | Ref. | Ø 4,5 mm | Ref. | Ø 5,0 mm | Ref. |
| Implant  | 3,5 x 7 mm | 8145 | 4,0 x 7 mm | 5904 | 4,5 x 7 mm | 6309 | 5,0 x 7 mm | 6315 | 4,5 x 7 mm | 5910 | 5,0 x 7 mm | 5917 |
| | 3,5 x 8,5 mm | 8146 | 4,0 x 8,5 mm | 5905 | 4,5 x 8,5 mm | 6310 | 5,0 x 8,5 mm | 6316 | 4,5 x 8,5 mm | 5911 | 5,0 x 8,5 mm | 5918 |
| | 3,5 x 10 mm | 5900 | 4,0 x 10 mm | 5906 | 4,5 x 10 mm | 6311 | 5,0 x 10 mm | 6317 | 4,5 x 10 mm | 5912 | 5,0 x 10 mm | 5919 |
| | 3,5 x 11,5 mm | 5901 | 4,0 x 11,5 mm | 5907 | 4,5 x 11,5 mm | 6312 | 5,0 x 11,5 mm | 6318 | 4,5 x 11,5 mm | 5913 | 5,0 x 11,5 mm | 5920 |
| | 3,5 x 13 mm | 5902 | 4,0 x 13 mm | 5908 | 4,5 x 13 mm | 6313 | - | - | 4,5 x 13 mm | 5914 | - | - |
| | 3,5 x 15 mm | 5903 | 4,0 x 15 mm | 5909 | 4,5 x 15 mm | 6314 | - | - | 4,5 x 15 mm | 5915 | - | - |
| Healing abutment  | 3,5 x 2 mm | 5922 | 3,5 x 2 mm | 5922 | 3,5 x 2 mm | 5922 | 3,5 x 2 mm | 5922 | 4 x 2 mm | 5925 | 4 x 2 mm | 5925 |
| | 3,5 x 3 mm | 5923 | 3,5 x 3 mm | 5923 | 3,5 x 3 mm | 5923 | 3,5 x 3 mm | 5923 | 4 x 3 mm | 5926 | 4 x 3 mm | 5926 |
| | 3,5 x 5 mm | 5924 | 3,5 x 5 mm | 5924 | 3,5 x 5 mm | 5924 | 3,5 x 5 mm | 5924 | 4 x 5 mm | 5927 | 4 x 5 mm | 5927 |
| Anatomic healing abutment  | 4 x 2 mm | 5928 | 4 x 2 mm | 5928 | 4 x 2 mm | 5928 | 4 x 2 mm | 5928 | 5 x 2 mm | 5934 | 5 x 2 mm | 5934 |
| | 4 x 3 mm | 5929 | 4 x 3 mm | 5929 | 4 x 3 mm | 5929 | 4 x 3 mm | 5929 | 5 x 3 mm | 5935 | 5 x 3 mm | 5935 |
| | 4 x 5 mm | 5930 | 4 x 5 mm | 5930 | 4 x 5 mm | 5930 | 4 x 5 mm | 5930 | 5 x 5 mm | 5936 | 5 x 5 mm | 5936 |
| | 5 x 2 mm | 5931 | 5 x 2 mm | 5931 | 5 x 2 mm | 5931 | 5 x 2 mm | 5931 | 6 x 2 mm | 5937 | 6 x 2 mm | 5937 |
| | 5 x 3 mm | 5932 | 5 x 3 mm | 5932 | 5 x 3 mm | 5932 | 5 x 3 mm | 5932 | 6 x 3 mm | 5938 | 6 x 3 mm | 5938 |
| | 5 x 5 mm | 5933 | 5 x 5 mm | 5933 | 5 x 5 mm | 5933 | 5 x 5 mm | 5933 | 6 x 5 mm | 5939 | 6 x 5 mm | 5939 |

ACCESSORIES

| | Ref. |
|--|------------------------------|
| Sterilization cassette Ocean | 1763 |
| MiniBox | 0526 |
| Torque wrench | 0295 |
| Screwdrivers | |
|  Screwdriver ISO 1797 S (048) Screwdriver ISO 1797 L (048) Screwdriver ISO 1797 XL (048) | 0644 0277 0645 |
|  Implant driver ISO 1797 EC/IC S (2.5) Implant driver ISO 1797 EC/IC L (2.5) Implant driver ISO 1797 W&H EC/IC (2.5) | 2693 0300 2692 |
|  Screwdriver gold screw ISO 1797 Screwdriver Rhein 83 abutment | 0263 |
|  Screwdriver grooved screw ISO 1797 | 0267 |
|  Screwdriver impression coping closed tray ISO 1797 | 0723 |
|  Screwdriver transepithelial abutment ISO 1797 S Screwdriver transepithelial abutment ISO 1797 L | 0328 0726 |
|  Screwdriver transepithelial angled abutment ISO 1797 S Screwdriver transepithelial angled abutment ISO 1797 L | 0804 0648 |
|  Implant driver ISO 1797 CC 3,5 S Implant driver ISO 1797 CC 3,5 L Implant driver ISO 1797 CC 4,1 S Implant driver ISO 1797 CC 4,1 L | 2981 2984 2982 2985 |


| | | | | |
|-------------------------------|---|------------------------|--|------------------------------|
| Screwdriver handle |  | Handle ISO 1797 S | 0791 | |
| |  | Handle ISO 1797 L | 0790 | |
| Implant handle |  | Implant handle S | 1878 | |
| |  | Implant handle L | 2891 | |
| Monoblock screwdrivers |  | Screwdriver S (048) | 0274 | |
| |  | Screwdriver L (048) | 0275 | |
| |  | Implant driver S (2.5) | 0278 | |
| |  | Implant driver L (2.5) | 0299 | |
| Drills |  | Pointed drill | 0188 | |
| |  | Twist drill | 1,6 x 7 - 15 mm | 2046 |
| |  | Pilot drill | 1,6 - 2,4 mm | 2047 |
| |  | Drill | 2,0 - 3,3 x 7 - 15 mm 2,2 - 3,8 x 7 - 15 mm 2,8 - 4,3 x 7 - 15 mm 3,2 - 4,8 x 7 - 13 mm | 2048 2049 2050 2051 |
| |  | Hard bone drill | 2,4 - 3,3 x 7 - 15 mm 3,0 - 3,7 x 7 - 15 mm 3,8 - 4,3 x 7 - 15 mm 4,1 - 4,7 x 7 - 13 mm | 3311 3312 3313 3314 |
| |  | Twist drill | 2,0 - 3,3 x 7 - 15 mm 2,2 - 3,8 x 7 - 15 mm 2,8 - 4,3 x 7 - 15 mm 3,2 - 4,8 x 7 - 13 mm | 8032 8033 8034 8035 |
| |  | Hard bone twist drill | 2,4 - 3,3 x 7 - 15 mm 3,0 - 3,7 x 7 - 15 mm 3,8 - 4,3 x 7 - 15 mm 4,1 - 4,7 x 7 - 13 mm | 8037 8038 8039 8040 |
| Screw tap |  | Screw tap | 3,5 mm 4,0 mm 4,5 mm 5,0 mm | 2687 2688 2689 2690 |

| | | | |
|---------------------------------|---|--------------|-----------|
| Direction indicator |  | 1,5 - 2,3 mm | 1810 |
| Drill extension ISO 1797 |  | | KI589B204 |

| | | | |
|---|---|----------------|------|
| Drill stops |  | Ø1,3-2,4 L7 | 6984 |
| | | Ø1,3-2,4 L8,5 | 6985 |
| | | Ø1,3-2,4 L10 | 6986 |
| | | Ø1,3-2,4 L11,5 | 6987 |
| | | Ø1,3-2,4 L13 | 6988 |
| | | Ø1,3-2,4 L15 | 6989 |
| |  | Ø2,8-3,3 L7 | 6990 |
| | | Ø2,8-3,3 L8,5 | 6991 |
| | | Ø2,8-3,3 L10 | 6992 |
| | | Ø2,8-3,3 L11,5 | 6993 |
| | | Ø2,8-3,3 L13 | 6994 |
| | | Ø2,8-3,3 L15 | 6995 |
| |  | Ø3,2-3,8 L7 | 6996 |
| | | Ø3,2-3,8 L8,5 | 6997 |
| | | Ø3,2-3,8 L10 | 6998 |
| | | Ø3,2-3,8 L11,5 | 6999 |
| | | Ø3,2-3,8 L13 | 7000 |
| | | Ø3,2-3,8 L15 | 7001 |
| |  | Ø3,6-4,3 L7 | 7002 |
| | | Ø3,6-4,3 L8,5 | 7003 |
| | | Ø3,6-4,3 L10 | 7004 |
| | | Ø3,6-4,3 L11,5 | 7005 |
| | | Ø3,6-4,3 L13 | 7006 |
| | | Ø3,6-4,3 L15 | 7007 |
|  | Ø4,2-4,8 L7 | 7008 | |
| | Ø4,2-4,8 L8,5 | 7009 | |
| | Ø4,2-4,8 L10 | 7010 | |
| | Ø4,2-4,8 L11,5 | 7011 | |
| | Ø4,2-4,8 L13 | 7012 | |
|  | Ø5,7 L7 | 7013 | |
| | Ø5,7 L8,5 | 7014 | |
| | Ø5,7 L10 | 7015 | |
| | Ø5,7 L11,5 | 7016 | |
| Drill Stops Kit for implants up to Ø4,5 |  | | 7959 |
| Drill Stops from Ø4,8 implants to Ø6,0 |  | | 7960 |


Instruments and surgical kit

Avinent supplies high-precision instruments designed to match implants and suitable for use with internal, external and conical connections. Avinent's two surgical boxes allow practitioners to choose the working system that best meets their needs.



The sterilization cassette OCEAN/ICEBERG stands out for its clear layout and attractive design, making it easy to use. The drilling sequence is clearly indicated by means of a simple colour code according to the diameter of the selected implant. All the pieces are placed in a sterilizable tray with a seethrough lid, giving a clear view of the interior.

The MiniBox is versatile, as it allows practitioners to select a specific sequence and take everything required for the surgical procedure with them in a small container. The box is sterilizable and can hold all the items needed for inserting prostheses.



Advisable torque

| TYPE | VALUE | |
|------------|-----------------------|---|
| Mechanical | 35 Ncm | Screw for single/Multiple abutment* |
| | 30 Ncm | Screw for angled titanium base** |
| | 30 Ncm | Screw for angulation correction G2 (max. 30°)** |
| | 20 Ncm | Screw for angulation correction G1 (max. 20°)** |
| | 15 Ncm | Screw transepithelial angled abutment |
| | 15 Ncm | Prosthetic screw for transepithelial |
| | 35 Ncm | Transepithelial abutment (Uniblock / 2 parts) |
| | 30 Ncm | LOCATOR® / LOCATOR R-Tx® Abutment |
| | 25 Ncm | RHEIN83® Abutment (OT Equator) |
| 15 Ncm | Temporary Abutment Ti | |
| Manual | ≈8-15 Ncm | Scan Abutment |
| | | Impression coping open/closed tray engaging |
| | | Healing Abutment Ti |
| | | Healing cap |

* Regarding: Titanium base, Cemented abutment and Cemented angled abutment, Castable CoCr Base. Included: Gold screw.

** Except Transep. 4.8 of M1.40 which is 15 Ncm.

Drilling speed biological protocol

| | rpm |
|------------------------|-------------|
| Guide drill | 800 - 1.200 |
| Drill ø 1,6 mm | 800 - 1.200 |
| Pilot drill | 600 - 800 |
| * Drill ø 2,0 - 3,3 mm | 150 - 300 |
| * Drill ø 2,2 - 3,8 mm | 150 - 300 |
| * Drill ø 2,8 - 4,3 mm | 150 - 300 |
| * Drill ø 3,2 - 4,8 mm | 150 - 300 |
| * Drill ø 2,4 - 3,3 mm | 150 - 300 |
| * Drill ø 3,0 - 3,7 mm | 150 - 300 |
| * Drill ø 3,8 - 4,3 mm | 150 - 300 |
| * Drill ø 4,1 - 4,7 mm | 150 - 300 |
| Screw Tap | 20 |

* The biological drilling at low revolutions (between 50 and 100 rpm), allows to collect autologous bone as shown in the scientific literature.

Maximum recommended torque for implant insertion: 45-50 Ncm

Maximum recommended speed for implant insertion: 20 rpm

Drill-bit length and marking


Avinent drills carry laser markings to improve visibility during osteotomy and follow a colour code according to the diameter of the implant. The marking corresponds to the length of the implant in crestal placement, but the distances are not absolute from the tip of the instrument to the mark. The length of the drill tip is not included in the depth mark, so this distance must be taken into account when planning treatment and in carrying out the osteotomy.

Drilling speed standard protocol

| | rpm |
|----------------------|-------------|
| Guide drill | 800 - 1.200 |
| Drill ø 1,6 mm | 800 - 1.200 |
| Pilot drill | 600 - 800 |
| Drill ø 2,0 - 3,3 mm | 200 - 400 |
| Drill ø 2,2 - 3,8 mm | 200 - 400 |
| Drill ø 2,8 - 4,3 mm | 200 - 400 |
| Drill ø 3,2 - 4,8 mm | 200 - 400 |
| Drill ø 2,4 - 3,3 mm | 200 - 400 |
| Drill ø 3,0 - 3,7 mm | 200 - 400 |
| Drill ø 3,8 - 4,3 mm | 200 - 400 |
| Drill ø 4,1 - 4,7 mm | 200 - 400 |
| Screw Tap | 20 |

Maximum recommended torque for implant insertion: 45-50 Ncm

Maximum recommended speed for implant insertion: 20 rpm



Accessories CC/CC.I

Implant handle


| | REF. |
|---|------|
| L | 2891 |
| S | 1878 |

Implant driver

| | REF. |
|-------|------|
| 3,5 S | 2981 |
| 3,5 L | 2984 |
| 4,1 S | 2982 |
| 4,1 L | 2985 |

Platform indicators

To make life easier for our costumers, Avinent screws attachments follow the color code of the implant platform.



* For any additional information and instrument maintenance instructions, please check www.avinent.com

Avent Biological drilling protocol

The Avent surgical drilling protocol for the ICEBERG system is suitable for all bone types. The system offers conical drills that adapt to the design of the implant. The finish of the drills makes it easier to locate the marks that indicate length during surgery. The biological drilling at low revolutions (between 50 and 100 rpm), allows to collect autologous bone as shown in the scientific literature.

The implant is designed so that the BAS surface treated part must be placed juxtaosseously and the polished part of 1.8 mm, transmucosally according to the clinical criteria.

These drills incorporate an optimized cut that guarantees better precision and efficiency in surgical procedures. In addition, DLC Treatment (Diamond-Like Carbon) technology has been integrated, which results in a type of coating that reduces wear, minimizes friction and extends the lifetime of the drills.

* Hard bone drills are identified with two color indicators.

Implant ø 3,5 mm

| | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,4 - 3,3 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 2048 | 3311 | 2687 | |
| | BONE TYPE III-IV | | | Option 1 | | Option 2 | |
| | BONE TYPE I-II | | | | | | |

Implant ø 4,5 mm

| | | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø2,8 - 4,3 mm | Drill ø3,8 - 4,3 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 2048 | 2049 | 2050 | 3313 | 2689 | |
| | BONE TYPE III-IV | | | | | | Option 1 | Option 2 | |
| | BONE TYPE I-II | | | | | | | | |

Implant ø 4,0 mm

| | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø3,0 - 3,7 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 2048 | 2049 | 3312 | 2688 | |
| | BONE TYPE III-IV | | | Option 1 | | Option 2 | | |
| | BONE TYPE I-II | | | | | | | |

Implant ø 5,0 mm

| | | | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø2,8 - 4,3 mm | Drill ø3,2 - 4,8 mm | Drill ø4,1 - 4,7 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 2048 | 2049 | 2050 | 2051 | 3314 | 2690 | |
| | BONE TYPE III-IV | | | | | | Option 1 | Option 2 | | |
| | BONE TYPE I-II | | | | | | | | | |

Aivent Standard drilling protocol

The Aivent drilling protocol for inserting ICEBERG system implants is based on the anatomical and morphological characteristics of the bone to be worked on, the aim being to achieve satisfactory primary stability and a good bone-implant contact percentage, key factors for good osseointegration. The diameter of the helicoidal drills in the system is related to the core of the implant to be inserted. Hence it is advisable to use the system's drilling sequence and instruments. The finish of the drills makes it easier to locate the marks that indicate length during surgery.

The implant is designed so that the treated part with the BAS surface being placed juxtaosseous and the polished 1.8 mm part transmucosally according to the clinical criteria.

These drills incorporate an optimized cut that guarantees better precision and efficiency in surgical procedures. In addition, DLC Treatment (Diamond-Like Carbon) technology has been integrated, which results in a type of coating that reduces wear, minimizes friction and extends the lifetime of the drills.

* Hard bone drills are identified with two color indicators.

Implant ø 3,5 mm

| | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,4 - 3,3 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 8032 | 8037 | 2687 | |
| | BONE TYPE III-IV | | | Option 1 | | Option 2 | |
| | BONE TYPE I-II | | | | | | |

Implant ø 4,5 mm

| | | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø2,8 - 4,3 mm | Drill ø3,8 - 4,3 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 8032 | 8033 | 8034 | 8039 | 2689 | |
| | BONE TYPE III-IV | | | | | | Option 1 | Option 2 | |
| | BONE TYPE I-II | | | | | | | | |


Implant ø 4,0 mm

| | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø3,0 - 3,7 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 8032 | 8033 | 8038 | 2688 | |
| | BONE TYPE III-IV | | | | Option 1 | | Option 2 | |
| | BONE TYPE I-II | | | | | | | |


Implant ø 5,0 mm

| | | | | | | | | | | |
|------|------------------|---------------|---------------------------|---------------------|---------------------|---------------------|---------------------|---------------------|-----------|---------|
| | Guide drill | Drill ø1,6 mm | Pilot drill ø1,6 - 2,4 mm | Drill ø2,0 - 3,3 mm | Drill ø2,2 - 3,8 mm | Drill ø2,8 - 4,3 mm | Drill ø3,2 - 4,8 mm | Drill ø4,1 - 4,7 mm | Screw tap | Implant |
| Ref. | 0188 | 2046 | 2047 | 8032 | 8033 | 8034 | 8035 | 8040 | 2690 | |
| | BONE TYPE III-IV | | | | | | Option 1 | Option 2 | | |
| | BONE TYPE I-II | | | | | | | | | |

Sterilization and packaging



Avinent implants are subjected to a sterilization process in accordance with the requirements of the CE mark for medical device. Avinent supplies its implants in packaging that is easily identified by a simple colour code. The vial cap comes in different colours according to the diameter of the implant and has a sticker indicating the type of connection and the diameter and length of the implant. The Avinent implant system is sold in a sterile blister pack, which ensures that the implant is fully protected until the given expiry date, so long as it is stored in suitable conditions. Avinent supplies labeling with all its products that ensure every item can be properly traced by means of stickers that can be used in the patient's clinical history or in any other documentation required.



Blister pack (front and back) Vial Opening the vial Implant Cover screw

Opening the packaging


- Open the box and take out the blister pack.
- Open the blister pack and leave the vial containing the implant in the sterile working area, avoiding contact with any item that is not sterile. The sterile blister pack must not be opened prior to use.
- Open the vial by gripping the coloured cap with one hand while holding the transparent part in the other.
- Fit the implant driver into the implant.
- Press to make a tight fit and raise the implant set slightly.
- The cover screw is inside the coloured cap.

Information on the box

Implant box label




Labels that make the product traceable



| | |
|---|---|
| Implant C.C.I. | ICEBERG conical connection indicator |
| ø 5.0 x 13 (4.1) | Diameter x length (platform) |
| Ti | Titanium |
| REF 8080 | Reference number |
| LOT AB12345678 | Batch code |
| YYYY-MM-DD | Expiry date |
| STERILE R | Sterile. Sterilization method: radiation |
| ifu.avinent.com | Consult instructions for use |
| Manufacturer | Manufacturer |
| CE 0197 | CE Mark NB num 0197. TÜV Rheinland |
| Rx only | (For the USA only) CAUTION: Federal law restricts this device to sale by or on the order of a physician |
| (01) 08435452001046 (10) AB12345678 (17) YYMMDD | UDI (Unique Device Identifier Carrier) |
| ⊗ | Do not re-use |
| ⊗ | Do not sterilize |
| ⊗ | Do not use if package is damaged and consult instructions for use |
| MD | Medical Device |
| Qty:1 | Quantity of implants included in the packaging |
| ⏳ | Date of manufacture |
| ⊗ | Simple sterile barrier system with protective packaging inside |

Avinent Guided surgery



Avinent Guided Surgery


- Surgical Guide Box 3DIEMME®
- Implant planning software with Real Guide®

CT or CBCT Scanners. Available libraries for:

- NewTom
- Carestream
- (More coming soon)

Surgical motors

Avinent offers a wide range of surgical motors designed for many clinical applications, produced by leading brands in the industry. Intended to simplify daily practice and for multiple users, they can be adapted to suit each practitioner's needs by allowing the entire drilling sequence to be customised. Our surgical motors are ergonomically designed, easy to use and significantly help to reduce tiredness. They give implantology specialists a surgical unit for everyday use in safely carrying out oral procedures. All the surgical motors have been designed with the aim of minimising size and weight, to improve the balance of the instrument and motor in the practitioner's hand and thereby reduce tiredness during long treatments.



Radiographic template

The radiographic template is a guide to help practitioners precisely select the appropriate implant diameter and length.

To adapt it to the proportions of the x-ray, the dimensions of the implant are shown at a scale 1:1 and 1.25:1.

Product index

| REF. | DESCRIPTION | REF. | DESCRIPTION |
|-----------|---|------|--|
| 0188 | Pointed drill | 878 | 5900 Implant CC.I 3,5 x 10 (3,5) |
| 0263 | Screwdriver gold screw ISO 1797 | 6 | 5901 Implant CC.I 3,5 x 11,5 (3,5) |
| | Screwdriver Rhein 83 abutment | | 5602 Implant CC.I 3,5 x 13 (3,5) |
| 0267 | Screwdriver grooved screw ISO 1797 | 6 | 5903 Implant CC.I 3,5 x 15 (3,5) |
| 0274 | Screwdriver S (048) | 7 | 5904 Implant CC.I 4,0 x 7 (3,5) |
| 0275 | Screwdriver L (048) | 7 | 5905 Implant CC.I 4,0 x 8,5 (3,5) |
| 0277 | Screwdriver ISO 1797 L (048) | 6 | 5906 Implant CC.I 4,0 x 10 (3,5) |
| 0278 | Implant driver S (2.5) | 7 | 5907 Implant CC.I 4,0 x 11,5 (3,5) |
| 0295 | Torque wrench | 6 | 5908 Implant CC.I 4,0 x 13 (3,5) |
| 0299 | Implant driver L (2.5) | 7 | 5909 Implant CC.I 4,0 x 15 (3,5) |
| 0328 | Screwdriver transepithelial abutment ISO 1797 S | 6 | 5910 Implant CC.I 4,5 x 7 (4,1) |
| 0526 | MiniBox | 6 | 5911 Implant CC.I 4,5 x 8,5 (4,1) |
| 0644 | Screwdriver ISO 1797 S (048) | 6 | 5912 Implant CC.I 4,5 x 10 (4,1) |
| 0645 | Screwdriver ISO 1797 XL (048) | 6 | 5913 Implant CC.I 4,5 x 11,5 (4,1) |
| 0648 | Screwdriver transepithelial angled abutment ISO 1797 L (0,48) | 6 | 5914 Implant CC.I 4,5 x 13 (4,1) |
| 5915 | Implant CC.I 4,5 x 15 (4,1) | 5 | 5917 Implant CC.I 5,0 x 7 (4,1) |
| KI589B204 | Drill extension | 7 | 5918 Implant CC.I 5,0 x 8,5 (4,1) |
| 0723 | Screwdriver impression coping closed tray ISO 1797 | 6 | 5919 Implant CC.I 5,0 x 10 (4,1) |
| 0726 | Screwdriver transepithelial abutment ISO 1797 L | 6 | 5920 Implant CC.I 5,0 x 11,5 (4,1) |
| 0790 | Handle ISO 1797 L | 7 | 5922 Healing abutment CC.I 3,5 x 2 |
| 0791 | Handle ISO 1797 S | 7 | 5923 Healing abutment CC.I 3,5 x 3 |
| 0804 | Screwdriver transepithelial angled abutment ISO 1797 S (0,48) | 6 | 5924 Healing abutment CC.I 3,5 x 5 |
| 1763 | Sterilization cassette OCEAN | 6 | 5925 Healing abutment CC.I 4,1 x 4 x 2 |
| 1810 | Direction indicator ø 1,5 - ø 2,3 | 7 | 5926 Healing abutment CC.I 4,1 x 4 x 3 |
| 1878 | Implant handle S | 7 | 5927 Healing abutment CC.I 4,1 x 4 x 5 |
| 2046 | Twist drill ø 1,6 x 7 - 15 | 7 | 5928 Anatomic healing abutment CC.I 3,5 x 4 x 2 |
| 2047 | Pilot drill ø 1,6 - ø 2,4 | 7 | 5929 Anatomic healing abutment CC.I 3,5 x 4 x 3 |
| 2048 | Drill ø 2,0 - ø 3,3 x 7 - 15 | 7 | 5930 Anatomic healing abutment CC.I 3,5 x 4 x 5 |
| 2049 | Drill ø 2,2 - ø 3,8 x 7 - 15 | 7 | 5931 Anatomic healing abutment CC.I 3,5 x 5 x 2 |
| 2050 | Drill ø 2,8 - ø 4,3 x 7 - 15 | 7 | 5932 Anatomic healing abutment CC.I 3,5 x 5 x 3 |
| 2051 | Drill ø 3,2 - ø 4,8 x 7 - 13 | 7 | 5933 Anatomic healing abutment CC.I 3,5 x 5 x 5 |
| 2687 | Screw tap 3,5 mm | 7 | 5934 Anatomic healing abutment CC.I 4,1 x 5 x 2 |
| 2688 | Screw tap 4,0 mm | 7 | 5935 Anatomic healing abutment CC.I 4,1 x 5 x 3 |
| 2689 | Screw tap 4,5 mm | 7 | 5936 Anatomic healing abutment CC.I 4,1 x 5 x 5 |
| 2690 | Screw tap 5,0 mm | 7 | 5937 Anatomic healing abutment CC.I 4,1 x 6 x 2 |
| 2891 | Implant handle L | 7 | 5938 Anatomic healing abutment CC.I 4,1 x 6 x 3 |
| 2981 | Implant driver ISO 1797 CC 3,5 S | 6 | 5939 Anatomic healing abutment CC.I 4,1 x 6 x 5 |
| 2982 | Implant driver ISO 1797 CC 4,1 S | 6 | 6984 Drill Stop Ø1.3-2.4 L7 |
| 2984 | Implant driver ISO 1797 CC 3,5 L | 6 | 6985 Drill Stop Ø1.3-2.4 L8.5 |
| 2985 | Implant driver ISO 1797 CC 4,1 L | 6 | 6986 Drill Stop Ø1.3-2.4 L10 |
| 3311 | Hard bone drill ø 2,4 - ø 3,3 x 7 - 15 | 7 | 6987 Drill Stop Ø1.3-2.4 L11.5 |
| 3312 | Hard bone drill ø 3,0 - ø 3,7 x 7 - 15 | 7 | 6988 Drill Stop Ø1.3-2.4 L13 |
| 3313 | Hard bone drill ø 3,8 - ø 4,3 x 7 - 15 | 7 | 6989 Drill Stop Ø1.3-2.4 L15 |
| 3314 | Hard bone drill ø 4,1 - ø 4,7 x 7 - 13 | 7 | 6990 CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L7 |
| 5184 | Screwdriver for angulation correcting G2 24 mm | 7 | 6991 CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L8,5 |
| 5185 | Screwdriver for angulation correcting G2 32 mm | 7 | 6992 CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L10 |
| 8145 | Implant CC.I 3,5 x 7 (3,5) | 4 | 6993 CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L11,5 |
| 8146 | Implant CC.I 3,5 x 18,5 (3,5) | 4 | 6994 CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L13 |

*The availability of the products in this catalog may vary depending on the country where you are based. Please contact your Avinent distributor for further information.

| REF. | DESCRIPTION | |
|------|---|----|
| 6995 | CORAL/OCEAN/ICEBERG Drill Stop Ø2,8-3,3 L15 | 8 |
| 6996 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L7 | 8 |
| 6997 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L8,5 | 8 |
| 6998 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L10 | 8 |
| 6999 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L11,5 | 8 |
| 7000 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L13 | 8 |
| 7001 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,2-3,8 L15 | 8 |
| 7002 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L7 | 8 |
| 7003 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L8,5 | 8 |
| 7004 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L10 | 8 |
| 7005 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L11,5 | 8 |
| 7006 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L13 | 8 |
| 7007 | CORAL/OCEAN/ICEBERG Drill Stop Ø3,6-4,3 L15 | 8 |
| 7008 | CORAL/OCEAN/ICEBERG Drill Stop Ø4,2-4,8 L7 | 8 |
| 7009 | CORAL/OCEAN/ICEBERG Drill Stop Ø4,2-4,8 L8,5 | 8 |
| 7010 | CORAL/OCEAN/ICEBERG Drill Stop Ø4,2-4,8 L10 | 8 |
| 7011 | CORAL/OCEAN/ICEBERG Drill Stop Ø4,2-4,8 L11,5 | 8 |
| 7012 | CORAL/OCEAN/ICEBERG Drill Stop Ø4,2-4,8 L13 | 8 |
| 7013 | OCEAN Drill Stop Ø5,7 L7 | 8 |
| 7014 | OCEAN Drill Stop Ø5,7 L8,5 | 8 |
| 7015 | OCEAN Drill Stop Ø5,7 L10 | 8 |
| 7016 | OCEAN Drill Stop Ø5,7 L11,5 | 8 |
| 7959 | Drill Stop Kit for implants up to Ø4.5 | 8 |
| 7960 | Drill Stop Set for implants from Ø4.8 to Ø6.0 | 8 |
| 8032 | OCEAN/ICEBERG Twist Drill Ø2.0-3.3 L7-15 | 14 |
| 8033 | OCEAN/ICEBERG Twist Drill Ø2.2-3.8 L7-15 | 14 |
| 8034 | OCEAN/ICEBERG Twist Drill Ø2.8-4.3 L7-15 | 14 |
| 8035 | OCEAN/ICEBERG Twist Drill Ø3.2-4.8 L7-13 | 14 |
| 8037 | OCEAN/ICEBERG Twist Drill Ø2.4-3.3 L7-15 | 14 |
| 8038 | OCEAN/ICEBERG Twist Drill Ø3.0-3.7 L7-15 | 14 |
| 8039 | OCEAN/ICEBERG Twist Drill Ø3.8-4.3 L7-15 | 14 |
| 8040 | OCEAN/ICEBERG Twist Drill Ø4.1-4.7 L7-13 | 14 |





Carretera de Navarcles, 107
Pol. Industrial Santa Anna I
08251 Santpedor (Barcelona) - España
T. (+34) 93 827 34 65
www.avinent.com - avinent@avinent.com

AVINENT AUSTRALIA
australia@avinent.com

AVINENT BENELUX
benelux@avinent.com

AVINENT CANADA
canada@avinent.com

AVINENT COLOMBIA
colombia@avinent.com

AVINENT CROATIA
croatia@avinent.com

AVINENT FRANCE
france@avinent.com

AVINENT KUWAIT
kuwait@avinent.com

AVINENT LITHUANIA
lithuania@avinent.com

AVINENT POLAND
poland@avinent.com

AVINENT PORTUGAL
portugal@avinent.com

AVINENT ROMANIA
romania@avinent.com

AVINENT TAIWAN
taiwan@avinent.com

AVINENT UKRAINE
ukraine@avinent.com

AVINENT UNITED ARAB EMIRATES
uae@avinent.com

AVINENT UNITED KINGDOM
uk@avinent.com

AVINENT USA
usa@avinent.com

AVINENT VIETNAM
vietnam@avinent.com

AVINENT SAUDI ARABIA
ksa@avinent.com



Management
System
ISO 13485:2016
www.tuv.com

